ENGLISH

⚠Caution: U.S. Federal law restricts this device to sale only by or on the order of a licensed dentist or physician.

Device Description

expanders are made in a variety of dimensions to accommodate differing patient anatomy. *expanders* have a wedge shape that is offered in maximum buccolingual thicknesses of 1.6, 2.0, 2.4, 2.8, 3.2, and 3.6 mm with lengths ranging from 9 mm to 15 mm. All *expanders* are 6 mm wide and manufactured from titanium alloy. *expanders* are provided in a non-sterile condition and must be cleaned and sterilized prior to use.

expanders provide the clinical benefit of expanding maxillary or mandibular bone that is at least 9-16 mm high and 3 mm wide. The patient target group for these devices are any persons requiring bone expansion. These devices may be used in conjunction with bone grafting and/or dental implant placement.

Indications for Use

expanders are intended for bone expansion in the maxilla and mandible.

Contraindications

Do not use *expanders* in presence of following:

- Patients who suffer from any medical conditions that make surgery inadvisable or may be otherwise deemed to be a contraindication by the treating dental practitioner,
- Titanium allergy,
- An insufficient quantity of bone volume (height and width) at the site, or
- A degree of bone mineralization not suitable for expansion.

Warnings and Precautions

- This device is restricted for sale only to or on the order of a licensed dentist or physician.
- The safe and effective use of *expanders* may be achieved only if these devices are used by qualified surgeons with appropriate training in bone surgery and expansion procedures.
- The safe and effective use of *expanders* is guaranteed only if devices are used as supplied per the instructions provided by Rex Implants. Modifications to the devices or disregard for these instructions for use may result in serious patient injury or death.
- Improve surgical technique including a preoperative study that does not properly evaluate the bone volume and bone quality may lead to bone damage and/or bone loss.
- Poor surgical planning including failure to reconcile the actual dimensions of surgical instruments to the sizes in radiographic measurements may result in extending the surgical site beyond the intended depth. This may result in permanent tissue damage that could cause hemorrhaging and/or permanent facial numbness.
- Failure to select proper *expander* size may result in tissue damage including bone fracture.
- To prevent accidental device aspiration and/or swallowing, tether a line to each device before use.
- Permanent numbness may occur if *expanders* contact nervous tissue.
- Trauma due to bone overheating and damage to adjacent tissues during the osteotomy must be minimized.
- Do not twist *expanders*. Torsion may result in bone and/or device fracture.
- Always begin insertion of *expanders* with a light force, increasing force only as necessary. Excessive force may result in bone damage, device fracture, patient dizziness and/or patient nausea.
- Do not reprocess devices for more than 50 separate cycles. Over-use of device may harm a patient, reduce surgical effectiveness, and/or lead to the device failure.
- In the event of device failure, ensure that all fragments are removed from the patient. Do not reuse device.
- Follow current local regulations and current facility procedure for the safe disposal of devices; devices should be cleaned and sterilized before disposal.
- Providing instructions to the patient is essential to ensure proper healing. The patient must be made aware of the limitations of the procedure, importance of oral hygiene, and avoidance of contraindications, and risk of adverse effects resulting from bone expansion surgery.

Potential adverse events

Potential adverse events must be communicated to the patient prior to surgery. Potential adverse events related to the use of *expanders* may include: wound dehiscence requiring bone graft; jaw bone fracture; the perforation of the following: maxillary sinus, inferior border of the mandible, labial and lingual bone walls, alveolar canal, and gingiva; abscesses; fistulas; supuration; inflammation; persistent pain; sensitivity reduction; paresthesia; hyperplasia; excessive bone reduction requiring surgery; systemic infections; nerve lesions or other nerve damage; and vascular lesions or hemorrhaging, which at times may be serious especially in patients undergoing treatment with anticoagulants and/or antiaggregants.

Failure to comply with the instructions contained in this document may cause harm to the patient, including risk of serious injury or death. Any serious incident that occurs in relation to these devices should be reported to Rex Implants Inc. and the relevant governing agency such as the US FDA or competent authority of the Member State in which the user and/or patient is established.

How supplied

expanders are provided non-sterile and must be cleaned and sterilized prior to use.

Cleaning/Reprocessing

expanders must be thoroughly cleaned before use. Reusable devices must be cleaned between uses. The cleaning processes should be carried out immediately after use in order to prevent the contaminants from drying onto the devices. While it is recommended that the following validated steps are included in a reprocessing protocol, the end-user bears the ultimate responsibility for the cleanliness of the device. These instructions are not intended for devices not manufactured by Rex Implants such as piezoelectric cutting devices.

Manual Cleaning/Disinfecting/Reprocessing Method

- Rinse device in cold, potable tap water (<43°C; <109°F) to remove debris and prevent blood coagulation.
- Prepare a solution of enzyme detergent* and potable tap water at pH 7, according to the manufacturer's instructions.
- Place the device in a clean container. Add a sufficient amount of enzyme detergent solution in the container to cover it completely.
- Leave the device immersed for 10 minutes at 40°C (104°F). This will reduce the organic residues. When immersed in the enzyme detergent solution, gently brush any holes, parts, and small grooves on the device using a toothbrush with soft nylon bristles until all visible contamination is removed.
- Place the device in an ultrasonic bath containing enzymatic detergent solution at 40°C (104°F) for at least 10 minutes. This will reduce the organic compounds on the devices. After soaking in the ultrasonic bath, gently brush any holes or grooves on the device.
- Under running, warm, potable tap water and using a toothbrush with soft nylon bristles, clean the device thoroughly without damaging the surface. Do this until all visible traces of soil are removed. Carry out the final rinse with distilled water.
- Finish cleaning and inspect the devices under an appropriate light source, paying attention to details that might house soil (e.g., laser etchings, holes, slots) and if necessary, repeat the cleaning cycle.
- Inspect all devices for any signs of wear and tear. Do not use any device whose integrity is visibly compromised.
- Dry the device in preparation for sterilization.

*Procedure validated with All-in-One 4 enzyme detergent.

Automated Cleaning/Disinfecting/Reprocessing Method

1. Perform a pre-cleaning by following steps 1 - 6 of the Manual Cleaning/Disinfecting/Reprocessing Method.

2. Lay the pre-cleaned device in a metallic tray and place it in thermosinfecter**.

NOTE: Place the devices in the washing machine so that dead zones do not arise and the water can properly drain. Also, make sure that the devices are properly held in place in the washing basket and cannot move during the washing process, as shocks could damage them.

⚠WARNING: Avoid overloading. Overloading compromises cleaning effectiveness.

3. Set the following sequence and parameters applicable to the cleaning cycle:

- 1 min, Pre-wash with cold, potable water;
- 5 min, Wash with alkaline detergent** at 55°C ± 2° (131°F);
- 1 min, Neutralization with proper solution** (1/3 cold water, 2/3 warm water);
- 1 min, Rinse with potable water (1/3 cold water, 2/3 warm water).

4. Disinfection:

- 5-min, Thermosinfect at 93°C (200°F) with demineralized water following the national requirements about A0 values;

- The automated thermosinfection was not experimentally tested. According to ISO standard 15883-1, Table B.1 (4) the thermosinfection at a temperature of 90°C [200°F] for 5 min, results in an A0 value of 3000.

⚠WARNING: Once the cleaning cycle in thermosinfecter has been completed, devices retain the washing temperature for a long time. When removing devices from the thermosinfecter, take proper precautions to avoid personal harm.

5. Inspect the devices under an appropriate light source, paying attention to details that might house soil (e.g., laser etchings, holes, slots), and if necessary, repeat the cleaning/disinfection cycle.

6. Inspect all devices for any signs of wear and tear. Do not use any device whose integrity is visibly compromised.

7. Dry the device in preparation for sterilization.

**Procedure validated with Miele PG8535 washer/disinfecter. Programm Miele DES-VAR-TD. Alkaline detergent: neodisher® FA (0,2 % v/v). Neutralizer liquid: neodisher® Z (0,1 % v/v)

Drying

Before starting the sterilization cycle, make sure that the device is thoroughly dry, both externally and internally. For this purpose, blow compressed air both externally and into/through any holes; this will prevent the onset of stains or halcos on the device.

Sterilization

expanders are provided non-sterile and must be cleaned and sterilized prior to use. These instructions are not intended for devices not manufactured by Rex Implants. Prior to sterilization, wrap all devices individually using standard self-sealing sterilization pouches that must be large enough not to stretch the seal or tear the packaging. Drying is to be performed inside the steam sterilizer in the prevacuum cycle at 132°C (270°F).

All sterilization phases must be performed by the operator in compliance to ANSI/AAMI/ISO 17665-1:2006, EN ISO 556-1:2001, ANSI/AAMI ST79:2017, and all appropriate requirements of the relevant governing body such as the US FDA or EU Member State(s). Care must be taken not to exceed the maximum load of the autoclave when sterilizing more than one instrument in the same cycle. Sterilization must be carried out using a pre-vacuum autoclave only. Any other sterilization methods must be avoided due to potential incompatibility with the materials used to manufacture the devices. At the end of the sterilization cycle, let devices cool down completely prior to use.

Sterilization validation has shown the following recommendations for sterilization to be effective to a SAL of 10⁻⁶.

Method: Steam

Cycle: Prevacuum for three cycles

Temperature: 132°C (270°F; tolerance 0°C to +3°C)

Minimum Exposure Time: 4 minutes

Minimum Drying Time: 20 minutes

Surgical Technique Synopsis

Step 1 – Preoperative Planning: Perform presurgical planning per standard surgical practices and techniques. This must include, but is not limited to, evaluation of bone quantity and quality for suitability for expansion, determination of optimal osteotomy size and shape, and selection of desired *expander* length. *expanders* are designed exclusively for the expansion of thin residual crestal bone at least 9 mm high and 3 mm wide.

Step 2: Create an osteotomy with a PIEZOSURGERY® cutting instrument that is at least 0.35 mm wide. Osteotomy depth should be 1 mm deeper than the desired *expander* length. A PIEZOSURGERY® device by Mectron is strongly recommended for use in this protocol. Piezosurgery has been demonstrated to improve surgical control and reduce the risk of tissue damage.

Step 3: Place a 1.6 mm *expander* with desired length in osteotomy. Insert flat end of device and maintain a vertical device orientation.

Step 4: Expand bone by applying force until the mechanical stop reaches the bone.

Step 5: Remove *expander* by applying force without bending or twisting device. Place a finger on top of the *expander* to stabilize it during removal.

Step 6: If additional expansion is desired, place a thicker *expander* in the osteotomy and repeat steps 4 and 5.

IFU-05 Rev 0, issued September 30, 2019. Additional information is available at <http://www.reximplants.com>.

Reference number and symbol	Title of symbol	Description of symbol per Standard*
5.1.1	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.
5.1.2	Authorized representative in the European Community	Indicates the Authorized representative in the European Community
5.1.3	Date of manufacture	Indicates the date when the medical device was manufactured
5.1.5	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
5.1.6	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
5.4.4	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
5.2.8	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
5.4.3	Consult instructions for use	Indicates the need for the user to consult the instructions for use
RxOnly 21CFR801.109(b)(1)	Prescription use only	Caution: U. S. federal law restricts this device to sale only by or on the order of a licensed dentist or physician
*Unless otherwise noted, reference numbers (e.g., 5.1.1) and descriptions from ISO 15223-1:2016, Medical Devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements, FDA recognized standard # 5-118		

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